

DNA Counterstain used by PSC

BRM and GD - 2/27/2019

NOTE – This protocol is appropriate for most IF slides and is required if images are obtained with the Aperio ScanScope FL. Read the reference documents carefully for more information.

Supplies Needed:

Reagent	Stock Concentration	Supplier	Reference
Hoechst 33258*	10 mg / mL	Molecular Probes #H3569	http://tools.lifetechnologies.com/content/sfs/manuals/mp21486.pdf
Deionized Water (dH ₂ O)		Purified in-house	
ProLong Gold Antifade Mounting Medium		Life Technologies #P36930	http://tools.lifetechnologies.com/content/sfs/manuals/mp36930.pdf

Procedure:

1. Prepare dilutions of Hoechst 33258. The stock solution is 10 mg/mL. Make a working solution of 100 µg/mL (1:100 dilution) in dH₂O and store in amber or light-resistant Eppendorf tubes at 4°C. Stock solutions in water are stable for six months when refrigerated.
2. At the time of staining, dilute the working solution (again 1:100 in dH₂O) for a final concentration of 1 µg/mL.
3. Add enough of the working solution to cover the tissue, then incubate at room temperature for 15 minutes.
4. Rinse slides three times (3x) in dH₂O, 5 minutes per rinse, using a light-resistant container.
5. Remove any excess water from around the tissue before mounting the coverslip. A folded Kimwipe works well for this purpose.
6. Add 1-3 drops of pre-warmed ProLong Gold mounting medium to the tissue. The quantity depends on the size of the tissue section. Then apply the coverslip[†].
7. Incubate the slides in the dark for 24 hours at room temperature. NOTE – The refractive index increases with curing time.
8. Seal coverslip with Nail polish and let dry if prolonged storage is needed or thick (>5 µm) sections are being prepared.
9. View the slides as usual.

* Hoechst 33342 or DAPI can be substituted for Hoechst 33258 but may require a slightly different dilution. The laboratory should perform a titration experiment to determine the optimal concentration for representative specimens.

[†] TPL has found that coverslips often have a residue that appears as blue artifact under fluorescent microscopy. This residue can usually be removed by a quick rinse with deionized water on both sides of coverslip before mounting.